



XERACOLOUR RANGE

CODE	Description	CI Numbers	pH	Density (g/cm ³)	Viscosity (KU@25°C)
XA	Xeracoulour Bright Yellow Interior	P.Y. 74	7.00 – 8.50	1.42 – 1.52	90 – 95
XALF	Xeracoulour Light Fast Bright Yellow	P.Y. 184	7.00 – 8.50	1.67 – 1.77	80 – 90
XAP	Xeracoulour Permanent Yellow	P.Y. 74/ 97	7.00 – 8.50	1.42 – 1.52	90 – 95
XAT	Xeracoulour Medium Yellow	P.Y. 83	7.00 – 8.50	1.30 – 1.40	90 – 95
XB	Xeracoulour Black	P.Bk. 7	7.00 – 8.50	1.40 – 1.50	95 – 100
XBPA	Xeracoulour Yellow Paste	P.Y. 74/53	7.00 – 8.50	1.65 – 1.75	75 – 85
XBPD	Xeracoulour High Strength Phthalo Green	P.G. 7	7.00 – 8.50	1.25 – 1.35	65 – 75
XBPE	Xeracoulour High Strength Phthalo Blue	P.B. 15.1	7.00 – 8.50	1.08 – 1.18	60 – 70
XBPH	Xeracoulour Red Paste	P.R. 112	7.00 – 8.50	1.20 – 1.30	70 – 75
XBPM	Xeracoulour High Strength Maroon	P.R. 122	7.00 – 8.50	1.06 – 1.14	80 – 90
XBPU	Xeracoulour Mid Yellow Paste	P.Y.65/83	7.00 – 8.50	1.26 – 1.34	70 – 80
XC	Xeracoulour Yellow Oxide	P.Y. 42	7.00 – 8.50	1.83 – 1.88	95 – 100
XD	Xeracoulour Phthalo Green	P.G. 7	7.00 – 8.50	1.48 – 1.55	95 – 100
XE	Xeracoulour Phthalo Blue	P.B. 15.1	7.00 – 8.50	1.45 – 1.55	95 – 100
AXF	Xeracoulour Red Oxide	P.R. 101	7.00 – 8.50	2.00 – 2.10	90 – 100
XF	Xeracoulour Red Oxide	P.R. 101	7.00 – 8.50	2.00 – 2.10	90 – 100
XG	Xeracoulour Light Fast Bright Orange	P.O. 73	7.00 – 8.50	1.43 – 1.53	85 – 95
XH	Xeracoulour Bright Red	P.R. 112	7.00 – 8.50	1.50 – 1.55	95 – 100
XHP	Xeracoulour Light Fast Bright Red	P.R.254	7.00 – 8.50	1.45 – 1.52	95 – 100
XJ	Xeracoulour Violet	P.V.23	7.00 – 8.50	1.55 – 1.65	95 – 100
XK	Xeracoulour White	P.W.6	7.00 – 8.50	1.90 – 2.00	95 – 100
XL	Xeracoulour Raw Umber	P.Y.42, P.R.101, P.BK. 7	7.00 – 8.50	1.56 – 1.65	95 – 100
XM	Xeracoulour Magenta	P.R.122	7.00 – 8.50	1.40 – 1.50	95 – 100
XOO	Xeracoulour Orange Iron Oxide	P.Y. 42	7.00 – 8.50	1.60 – 1.70	85 – 95
XQ	Xeracoulour Yellow Light Fast	P.Y.184	7.00 – 8.50	1.55 – 1.65	85 – 95